

---

# TEAM 205 – ROBOTIC ARM TO ASSIST MOTORIZED CHAIR USERS

Industry Sponsor: PPMD, Mr. Keith Van Houten

Faculty Advisor: Prof. Jim Hess


Team Members: Brooke Harrington, Cavan Moriarty, Omar Abuljobain,  
Justin Lee, and Andrew Baum

---



# Agenda

---



# Problem Background – Context for Future Design

## Background

Duchenne Muscular Dystrophy (DMD) is a genetic disorder that mainly affects boys and causes progressive muscle weakness as they age

People with DMD require a motorized wheelchair for general movement and a caretaker or assistive device for other tasks throughout the day



[Figure 1: Permobil M3 Corpus](#)

## Problem

As people with DMD get older and their muscles are significantly weakened, they require more sophisticated assistive devices that aid upper body mobility

These devices are often robotic aids that function as arms and hands which require complicated controls. Additionally, current assistive robotics on the market are priced >\$60000 and are often not covered by insurance.



[Figure 2: Kinova JACO Robotic Arm. Highly complex arm that requires toggle between which part is being controlled, which includes the fingers](#)

## Value Proposition

An affordable, easy to control robotic aid allows more people with DMD to have access to a device that provides range of motion (ROM) assistance for daily tasks that include pushing elevator buttons or grabbing a drink to bring to the user



# Architecture - How the System Functions

## Mechanical Systems



### Mounting Structure



### Case



## Electrical/Software Systems



**Function:**

- Connects to the Permobil Unitrack for easy installation/removal
- Adjustment of system position for comfortability
- Mounting for the case

**Function:**

- Provides structure and mounting for the Waveshare robotic arm
- Contains and protects the electronics
- Mounting for control devices and interface (joystick, buttons, switches, and LEDS)

Control Panel (JSTK, BTN, SW, LED)

Analog and Digital Signals

ESP32 Board on Robot

JSON Document

Adafruit Feather

Motor Control

End Effector Coordinate Control

# Sys Obs – Two Key Functions Desired by End Users



## Pushing Elevator Buttons:

- System is capable of depressing common elevator buttons, both call and destination buttons
- A “ready” position can be saved so that a single button push will bring the end effector near the location of the buttons, then, the joystick will be used to fine tune the position and hit the correct button





## Automatic Drink Functions:

- Automatic drink function where the user programs the path the end effector should take to the drink and user. This function can be later recalled through the push of a button
- End effector control through the joystick buttons which allows the user to retrieve a drink in a specific portable cupholder

# Initial Concepts and Down Selection Process

Criteria for Screening and Scoring	Weight
Operating Radius	10%
Number of Inputs	15%
Number of Actuators (DOF)	10%
Size	15%
Cost	10%
Safety	15%
Failure Modes	10%
Integration into Chair	10%
Ease of Installation	5%




## Final Concept and Key Features

- Combination of some of the concepts developed by members
- Joystick, robot, battery and all accompanying electronics contained in the same box that is easily removed
- Simple robot, max five degrees of freedom
- Mounts near/under the armrest on the Permobil wheelchair and uses the Unitrack as the attachment point

# Key Factors, Engineering Analysis

## Structure:

- Initial FEA analysis completed for the vertical rails to ensure system would not fail under its own weight
- FEA showed a max stress of 48.436 MPa (FOS 5.67 when compared to yield strength of anodized aluminum)



## End Effector Accuracy:

- Calculations for accuracy at different points along the robotic arm based on advertised motor accuracy
- Calculations showed that the selected robot would be able to accurately push elevator buttons compliant with ADA standards





Figure 3: ADA standards for elevator buttons

# CAD Breakdown: Overall System

## Key Components


- Unitrack Mounting System
- Electronics Case
- Control Panel
- 3D Printed Armrest
- Waveshare Ro-Arm



## Key Features

- Permobil Wheelchair Mounting
- Horizontal and Vertical Adjustment
- Protected Electronics
- Intuitive and Familiar Controls
- Arm support
- Interactions with Surroundings

# CAD Breakdown: Case




## Key Components

Outer Acrylic Plates:


- Control Panel
  - Joystick
  - Switch (x2)
  - Buttons (x5)
  - LED (x3)
- Lithium Ion Battery and Battery Case
- Waterproofing Well

3D Printed Frame:

- Upper French Cleats
- Armrest
- Protoboard Assembly mounting
  - Adafruit Feather
- Waveshare Ro-Arm mounting
  - ESP 32



# CAD Breakdown: Mounting System




## Key Components

### Mounting Assembly:


- 1/4" Aluminum Mounting Plate
- 45° T-Slotted Rails
- Wing-Bolts
- Square Nuts

### Connecting Assembly:

- 1/8" Aluminum Connecting Plate
- Hex Head Bolts
- Slot for 3D Printed Key
- Lower French Cleats



# CAD Breakdown: Safety Features




## Key Components

### Mechanical Stop:

- 1/8" 3D Printed Stopping Plate
- Standoff


### Waterproofing:

- 1/8" 3D Printed Well
- 3D Printed Battery Cover
- 3D Printed End-Effector Cover



# Build Progress: Physical Systems, Key Features

Waveshare Robotic Arm with Waterproofing Solutions



Case Containing Electronics and User Interface



Mounting Structure




User Interface (Joystick, Buttons, Switches, LED)

Armrest

Vertical Adjustments Rails

# Build Progress: Physical Systems - Mounting




- The lower diagonal rail was cut to remedy safety concerns
- There was concern for the stability of the mounting schematic on the Jazzy wheelchair (vibrations) but proved to be stable on the Permobil through sponsor testing
- There is current concern for the reliability of the 3D-printed French cleat key - broken multiple times despite solid infill




# Build Progress: Physical Systems - Case





- Updated user interface with joystick shifted right
- 2 LEDs removed from user interface due to electronics complications
- Added armrest for additional support while in use
- Mounting frame with cutout for feather to slide into place under the ESP 32 headers on the robotic arm
- Nut inserts in frame for mounting acrylic plates and French cleats



# Build Progress: Electronics Systems




- The protoboard is the main connection point between the feather, control systems, and ESP32 on robot
- All control systems are soldered directly to the protoboard. Using separate connectors posed issues for long term robustness
- Pigtail cable that powers the robotic arm from the battery. Robotic arm powers the feather and all control systems



# Build Progress: Waterproofing Systems



- Added a “cap” to the end effector to shield from rain (clips into the LED mount)
- Sleeve to cover the “elbow” and shoulder joint, attaches to the carbon fiber tube that goes from the elbow to the end effector
- “Well” used under the stopping plate to keep water from falling onto the ESP 32 board under the robot



# Test Results

## Test 1: End Effector Max Length Accuracy

Addresses TR 1.3.2: The system shall have an accuracy of +/-10mm

### Test Equipment:

- Calipers, Tape measure, Whiteboard, Whiteboard marker, Tape, Wheelchair

### Summary of Test Procedure:

- Assemble test system onto wheelchair and attached marker onto end-effector with tape
- Mark 35" and 48" lines on whiteboard to replicate ADA standard button height dimensions
- Navigate wheelchair to the whiteboard so the user is perpendicular to the surface
- Navigate arm with marker and mark three points at 35" and 48" mark
  - This will be done at the lowest increment by slightly shifting the joystick past its dead zone
- Measure horizontal (X,Y) and vertical (Z) distances and compare to expected value

Expected Value: <10 mm (0.39")

Outcome Value: ~5-14 mm  
(0.19-0.55")

### Outcomes/Conclusions/Recommendations:

- Measured values fell below and above expected value
- The cause of fluctuation was the rate at which the joystick function is looped
- Meets the button standard of 0.75" set by ADA, code revision is still required



# Test Results

## Test 2: End Effector Max Load Capacity

Addresses TR 1.1.4: The arm should be able to bring a drink weighing 470g to the user

### Test Equipment:

- Drink (liquid, container, and portable cupholder), Scale, and User interface (drinking function)

### Summarized Test Procedure:

- Program drinking function (min. and max. range)
- Measure/record weight of drink: liquid, container, and portable cupholder
- Position drink at the initial position of the drink function, then initiate the drink function
- Observe/record the system behavior; determine if the system is stable or unstable
- Repeat steps and increase weight by 0.1kg until system is determined to be unstable

Expected Value: Min. and Max.  
(500g)

Actual Value: Min. (600g) and  
Max. (500g)

### Outcomes/Conclusions/Recommendations:

- This test was performed with clamp and table; recommend re-running test mounted to wheelchair. Similarly, test was performed without end effector cover
  - The robotic arm is capable of fulfilling TR.1.1.4



# Test Results

## Test 3: Restricted Motion Compliance

Addresses TR 5.2.0: The robot shall have physical and program safeguard to prevent it from reaching undesired locations

### Test Equipment:

- Robotic arm system, Monitoring laptop

### Summarized Test Procedure:


- Power on the system
- Initiate arm motion
- Test positive and negative arm limits (min. and max. range along arc)
- Test edge case use
- Repeat steps 3-6 at medium and maximum speed

Expected Value: Remained within +/- 30 degrees

Actual Value: Remained within +/- 30 degrees

### Outcomes/Conclusions/Recommendations:

- System passes test by not entering safe zone
- Physical motion restrictor was not utilized, digital motion restrictor sufficient for test



# Test Results

## Test 4: Installation Time

Addresses TR 3.0.0: The arm shall be installed in less than one hour and be operational within 30 minutes

### Test Equipment:

- Phone stopwatch, mounting system, electronics box

### Summarized Test Procedure:

- Start stopwatch
- Attach mounting system to Unitrack
- Mount electronics box to French cleats of mounting system
- Place key in slot
- Stop stopwatch

Expected Value:


30 min

Actual Value:

58.52 sec

### Outcomes/Conclusions/Recommendations:

- System passes the test
- Significant overprediction due to many design simplifications



# Test Results

## Test 5: Maintenance Time

Addresses TR 3.1.0: The system shall have maintenance work completed in less than three hours

### Test Equipment:

- Phone stopwatch

### Summarized Test Procedure:


- Begin timer, remove both side panels and the diagonal panel that contains the control interface
- Inspect the 3D printed frame and all fasteners. Inspect the wires of the control interface, and all solder joints/heat shrink/possible loose connection points
- Inspect the mounting structure visibly and by touch
- Reassemble the system
- Power on the system to ensure normal operation
- End timer, record time

Expected Value: 1 hr

Actual Value: 22 min 05 sec

### Conclusions/Recommendations:

- This maintenance inspection did not require any fixing of parts
- Recommend re-running this test after significant use or next maintenance period ~ 6 months



# Test Results – Meeting with Colin



Colin using the joystick to get a drink  
(Note this is not the automatic drink function – we were unable to test due to mounting interferences)



Colin going through door (door is larger than standard, but is ADA compliant)



Colin testing pushing buttons on lift outside – end effector sometimes slipped off the button



# Test Results –Colins Feedback



## Joystick Position:

- The initial design in TRR was developed under assumed dimensions of the chair, since each person's chair is different
- Fitment testing on Colins chair showed some initial assumptions were incorrect
- Joystick was moved right to align closer with the armrest on Colins chair

## Joystick Selection:

- Toggling the end effector between open and close was not ideal
- New joystick has two buttons on top to pitch the end effector to specific angles, other controls remained the same



## Case Changes:

- In addition to moving the joystick over, added an armrest for additional support
- Added padding on armrest for comfort
- Battery remains accessible through cutout underneath

# Status of Requirements – Updates since TRR

Technical Requirement ID	Requirement	Verification Method	Status
TR 1.1.0	The robot shall be controlled manually by a joystick and should have the ability to program automatic functions accessed through a button push	Test	Completed
TR 1.2.0	The robot shall have no more than six degrees of freedom to reduce the complexity of controls	Inspection	Met with robot arm choice
TR 1.3.0	The robot shall have an operating radius that spans a minimum of 5 inches to a maximum of 2 feet from mount location with 180 degrees of rotation total	System Test	Met with robot arm choice
TR 1.4.1	The battery that powers the robot arm shall last 24 hours between charges	Analysis	Completed

# Status of Requirements – Updates since TRR

---

Technical Requirement ID	Requirement	Verification Method	Status
TR 1.3.1	The device should be able to depress elevator buttons and the handicap door button between 35" and 48" as specified by ADA standards	Analysis	Completed
TR 1.3.2	The system shall have an accuracy of +-10mm while at full extension	Analysis and test	Analysis Completed – Testing Completed – Revision/Retest Needed
TR 3.0.0	The robot shall be installed in less than one hour and should be operational within thirty minutes of installation	Test	Completed
TR 5.3.0	The robot shall have safeguards against unexpected de-energization	Inspection	Completed

# Conclusions

## Meeting User Requirements:

- The team believes we have met most of the user requirements
- Need additional user testing and feedback to verify this assertion
- Colins feedback suggest the project is on the correct path: new user requirements can be added for project continuation

## Marketability:

### Team 205

~\$2,000

4 DoF

Detachable System

Open Source

### Market

~\$60,000

6 DoF

Fixed System

Built in Software

## Overall System Performance:

- Team is satisfied with performance and proud of the effort that went into developing a functional system
- Though satisfied, refinements can be made that would improve multiple facets of the system

# Lessons Learned

## What Went Well:

### Design Process:

- Constant feedback through many meetings
- Continuously evolving design
- End user's needs met

### Team Improvement:

- Significant refinement since PDR
- More information-quicker progress, more productivity

### Teamwork:

- Equal division of tasks
- Gain new skills

## Where We Fell Behind:

### After PDR:

- CAD Lacking Detail
- Missing Context
- Misunderstandings

### After CDR:

- Lacking electronics design
- Lacking a plan for interacting with the robot (microcontroller)
- Minimal code

### After TRR:

- Waterproofing

## Important Realizations:

### Design Process:

- No standard configuration – difficult design process
- Trying to overaccommodate leads to delays and confusion

### Constant Feedback:

- Helpful during design process
- Allowed for continuous refinement
- Led to user satisfaction

### Current State:

- Room for improvement
- Time constraints limit optimizations

# Future Work

## Future Work for our Team:

- Work on the report
- Finalize any last aspects of the project we feel need a little more work (very small tasks)
- Working on the final deliverables for Mr. Keith from PPMD
- GitHub – giving as much as we can so this can continue as an open-source project

## Things We Would Change Given More Time:

### Robotic arm choice:

- Initial analysis showed the Waveshare was a fair option for our project
- After working with the Waveshare, other options may be better, more research needs to be done

### Coding:

- Coding could use refinement and optimization
- None of the members have any formal education in robotics coding


## Recommendations for a Future Team:

- Maybe try a new arm – learn from our choice of arm and its capabilities to influence their own decisions
- Continue frequent discussion with the PPMD team – it was a huge help with the design process
- The case can be reworked, as well as the mounting mechanism, which can likely be strengthened and motorized for ease of positioning
- There's different, likely better ways to meet the customer needs that should be explored

# One Last Project Goal We Want to Cover

## Open Source:

- We are populating a GitHub with our CAD, code, and manuals files
- Updating a google drive with our tutorial videos, and other images from the duration of the project
- Developing a BOM (GitHub) for all materials in the final design, and estimates for manufacturing from third party companies
- QR codes for GitHub and images on expo poster



GitHub



Photos


# THANK YOU FOR LISTENING!

Before questions, the team would like to open the presentation to Mr. Keith or Colin via zoom to discuss more on DMD, PPMD and the long-term goals of this project




```
383 }
384 }
385
386 void programstow2() { //get water
387     // Start the flash memory read write
388     Preferences SAVEDPOSITIONS;
389     SAVEDPOSITIONS.begin("SAVEDPOSITIONS", false);
390
391     // Grab the current pos from JSON feedback and deserialization
392     StaticJsonDocument<200> initdoc;
393     initdoc["t"] = 105;
394     serializeJson(initdoc, Serial1);
395     Serial1.println();
396     delay(100);
397     String response = Serial.readStringUntil('\n');
398     StaticJsonDocument<200> responseDoc;
399     DeserializationError error = deserializeJson(responseDoc, response);
400     if (error) {
401         curxpos = responseDoc["x"];
402         curypos = responseDoc["y"];
403         curzpos = responseDoc["z"];
404         curtpos = responseDoc["t"];
405     }
406
407     //Update the state to Program Stow
408     bool running = true;
409
410     // While in ProgramStow, enable joystick control if no button is being pushed, and when pushed
411     // record the position values and write to EEPROM to be called later
412     while (running) {
413         if (digitalRead(D4) == LOW){
414             float stow2x = curxpos;
415             float stow2y = curypos;
416             float stow2z = curzpos;
417
418             digitalWrite(D0, HIGH);
419             SAVEDPOSITIONS.putString("Stow2XPosition", String(stow2x));
420             SAVEDPOSITIONS.putString("Stow2YPosition", String(stow2y));
421             SAVEDPOSITIONS.putString("Stow2ZPosition", String(stow2z));
422
423             delay(500);
424             digitalWrite(D0 ,LOW);
425
426             SAVEDPOSITIONS.end();
427
428             running = false;
429         }
430         else {
431             jstkcntl();
432             ledblinker();
433         }
434     }
435 }
```

# Appendix – Assembly Manual

Step 33: Replace the button head and locknut that fits into the rails with an M6 x 0.75, 12mm long bolt and square nut.



Step 34: Slide the case's French Cleats onto the mounting assembly's French Cleats.



Enjoy!

# Blooper Slide

---



# Bill of Materials

Team Number	205
Team Name	Robotic Arm to Assist Motorized Chair Users
Total Cost	\$ 1,068.92

Assist Motorized Chair Users

## Bill of Materials

Taxes and Shipping NOT included

Level-1: Electronics  
 Level-2: Waterproofing  
 Level-3: Mounting  
 Level-4: 3D Printed  
 Level-5: Hardware  
 Level-6: Other

BOM Level	Part Name	Part Number	Part Details	Unit Cost	Quantity (Ea./Pkg.)	Total Cost	Procurement Details	Additional Notes
Level-1	Waveshare Roarm	RoArm-M2-S		\$ 223.99	1 Ea.	\$ 223.99	<a href="#">Amazon</a>	Direct Order
Level-1	Joystick	HE2 SERIES 3 AXIS	Handle(two buttons), Limiter(square), Outout(0)	\$ 290.07	1 Ea.	\$ 290.07	<a href="#">Ruffy Controls</a>	Direct Order
Level-1	Buttons	BOC2G4RC4Q	Size(momentary) and Color(silver)	\$ 12.49	5 Pkg.	\$ 62.45	<a href="#">Amazon</a>	Direct Order
Level-1	Switches	BOCKXRHJPF	Type(5pcs switch + switch holder)	\$ 8.09	2 Pkg.	\$ 16.18	<a href="#">Amazon</a>	Direct Order
Level-1	Wires	1311	22AWG Solid Core	\$ 15.95	1 Ea.	\$ 15.95	<a href="#">Adafruit</a>	Direct Order
Level-1	Heat Shrink Tubing	asd-123	Heat Shrink Tube Kit with Box(5 colors/12 Sizes)	\$ 6.99	1 Pkg.	\$ 6.99	<a href="#">Amazon</a>	Direct Order
Level-1	Male-Male Plug Cable	B0BN2YX895	Voltage(12V) and Input Current(1.5A)	\$ 8.99	1 Pkg.	\$ 8.99	<a href="#">Amazon</a>	Direct Order
Level-1	Battery	PB120B1	12V Lithium ion	\$ 8.99	1 Pkg.	\$ 8.99	<a href="#">Amazon</a>	Direct Order
Level-2	Sleeve Cover	B0D6Y2R9LZ	Color(black+blue) and Material(PU leather)	\$ 9.56	1 Pkg.	\$ 9.56	<a href="#">Amazon</a>	Direct Order
Level-2	Hole Plugs	suiwotin - 220105 - 7	Color(white) and Size(1/4")	\$ 6.99	1 Pkg.	\$ 6.99	<a href="#">Amazon</a>	Direct Order
Level-2	Silicon Gasket Maker	BOCXPKXLCCH	Size(3 Fl. Oz.) and Style (Ultra Black)	\$ 14.96	1 Pkgs.	\$ 14.96	<a href="#">Amazon</a>	Direct Order
Level-3	Sliding Plate		Thickness(1/4") and Size(1x1ft)	\$ 52.42	1 Ea.	\$ 52.42		Additional Manufacturing
Level-3	Connecting Plate		Thickness(1/8") and Size(1x1\2ft)	\$ 26.47	1 Ea.	\$ 26.47		Additional Manufacturing
Level-3	Fixturing T-Track	1850A11	Length(18")	\$ 9.83	2 Ea.	\$ 19.66	<a href="#">McMaster</a>	Direct Order
Level-3	French Cleats	2152A11	Material(Al.)	\$ 8.24	1 Pkg.	\$ 8.24	<a href="#">McMaster</a>	Direct Order
Level-4	3D-Prints		Material(ABS)	\$ 95.56	1 Pkg.	\$ 95.56		Additional Manufacturing
Level-5	Heat-Set Inserts for Plastic	94459A340	Brass, 10-32 Thread Size, 1/4" Installed Length	\$ 12.36	1 Pkg.	\$ 12.36	<a href="#">McMaster</a>	Direct Order
Level-5	316 Stainless Steel Button Head Hex Drive Screw	98164A178	Corrosion-Resistant, 10-32 Thread Size, 1/2"	\$ 9.82	1 Pkg.	\$ 9.82	<a href="#">McMaster</a>	Direct Order
Level-5	Male-Female Threaded Hex Standoff	91075A467	ainless Steel, 3/16" Hex, 5/16" Long, 6-32 to 6-32	\$ 3.79	1 Ea.	\$ 3.79	<a href="#">McMaster</a>	Direct Order
Level-5	18-8 Stainless Steel Pan Head Slotted Screws	91792A144	6-32 Thread Size, 1/4" Long	\$ 5.51	1 Pkg.	\$ 5.51	<a href="#">McMaster</a>	Direct Order
Level-5	Lock Nut for T-Track	1850A26	Steel, 3/4"	\$ 3.38	2 Ea.	\$ 6.76	<a href="#">McMaster</a>	Direct Order
Level-5	Low-Strength Steel Square Nut	96887A319	Class 5, Zinc-Plated, M6 x 1 mm Thread	\$ 14.40	1 Pkg.	\$ 14.40	<a href="#">McMaster</a>	Direct Order
Level-5	Alloy Steel Socket Head Screw	91290A017	Black-Oxide, M2 x 0.4 mm Thread, 10 mm Long	\$ 15.56	1 Pkg.	\$ 15.56	<a href="#">McMaster</a>	Direct Order
Level-5	High-Strength Steel Nylon-Insert Locknut	94645A102	ss 10, Zinc Plated, M5 x 0.8 mm Thread, 5 mm H	\$ 16.62	1 Pkg.	\$ 16.62	<a href="#">McMaster</a>	Direct Order
Level-5	Low-Strength Steel Nylon-Insert Locknut	90633A007	Thin-Profile, Zinc-Plated, 6-32 Thread Size	\$ 10.20	1 Pkg.	\$ 10.20	<a href="#">McMaster</a>	Direct Order
Level-5	Stainless Steel Wing-Head Thumb Screw	92625A413	M6 x 1mm Thread Size, 12mm Long	\$ 4.44	4 Ea.	\$ 17.76	<a href="#">McMaster</a>	Direct Order
Level-6	Acrylic Panel		hite Opaque), Dimension(12"x12"), and Thickne	\$ 76.68	1 Pkg.	\$ 76.68		Additional Manufacturing
Level-6	Portable Cupholder	BOCKGM1G9M	Color(Black Reflective)	\$ 11.99	1 Ea.	\$ 11.99	<a href="#">Amazon</a>	Direct Order
<b>Total</b>							\$ 1,068.92	

